Biotest MMB Prep Campylobacter

Manual

Version 06.09

Biotest MMB Prep Campylobacter (2x 50 Preps)
Art. No. 944 511

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hycon@biotest.de

Introduction

Campylobacter jejuni, lari, coli are the most frequent causes of bacterial colon infections in humans worldwide. The transmission of these bacteria occurs mainly by undercooked food or by recontamination of food products, especially by poultry and ready-to-eat food. The most important *Campylobacter* species is the thermo-tolerant *C. jejuni*. Additionally, relevant are *C. coli* and *C. lari*, which are also associated with affection in human colon.

The real-time PCR method, in contrast to classic microbiological culture methods, enables the rapid and sensitive testing of raw material and processed food for absence of *Campylobacter*.

The **Biotest MMB Prep Campylobacter** kit can be used for template preparation for PCR analysis.

Intended Use

The kit is intended to be used for the isolation of *Campylobacter* DNA from enrichment cultures of food samples..

Kit Components (50 preparations) and Storage

- 1x Lysis Buffer (20 ml)
- 1x Binding Buffer (10 ml)
- 1x Pre-Wash Buffer (60 ml, after adding 30 ml ethanol (absolute, not included in the kit))
- 1x Wash Buffer (60 ml, after adding 42 ml ethanol (absolute, not included in the kit))
- 1x Elution Buffer (10 ml)

- 2x Receiver Tube 2 ml (50 x)
- 1x Receiver Tube 1.5 ml (50 x)
- 2x Spin Filter (50 x)

All kit components are stored dry and at room temperature (14-25 °C).

Additional required reagents

Ethanol (absolute, >96 %)

Additional required equipment and disposables

- Microcentrifuge
- Thermomixer / heating block (100 °C)
- Adjustable pipettes
- Sterile DNA-free filter tips
- 1.5 ml and 2 ml microcentrifuge tubes

Quality Control

According to the Quality Management System, each lot of Biotest MMB Prep Kits is tested against predetermined specifications to ensure consistent product quality.

Warning and Precautions

Not for *in vitro* diagnostic use.

- Caution:** The test must be performed by well-trained and authorized laboratory technicians.
- Caution:** All reagents should be handled in accordance with good laboratory practice taking appropriate precautions.
- Caution:** Biological samples can contain infectious agents. Follow all appropriate precautions, e.g., the regular check of gloves for contamination with sample material.
- Caution:** Do not use reagents past the expiration date printed on the label.
- Caution:** Do not use reagents with any evidence of turbidity or microbial contamination.

For more information, please consult the appropriate material safety data sheets (MSDS). These are available online in PDF format at www.biotest.de. In the section Microbiological Tests & Systems you can find, view, and print the MSDS for each Biotest MMB Prep kit and kit component.

Protocol

Before using the kit for the first time add 30 ml ethanol absolute to the Pre-Wash Buffer and mix thoroughly.

Before using the kit for the first time add 42 ml ethanol absolute to the Wash Buffer and mix thoroughly.

Before each isolation procedure the required amount of Elution Buffer is transferred to a reaction tube and pre-warmed to 65 °C.

Preparation of starting material

1. Transfer 1 ml of the enrichment culture under sterile conditions into a 1.5 ml reaction tube (not provided with the kit).
2. Centrifuge for 5 min at 14,000 x g.
3. Remove and discard supernatant by pipetting.

Lysis of starting material

4. Add 400 µl* of **Lysis Buffer**, mix briefly and close the reaction tube tightly.
5. Incubate under continuous shaking for 10 min at 99 °C.

* The volume depends on the specific sample matrix. It is possible to vary the volume up to 1.0 ml. For further questions please send contact to hycon@biotest.de.

Setting of optimal binding conditions

6. Centrifuge for 1 min at 14,000 x g.
7. Place a **Spin Filter** into yellow **2 ml Receiver Tube**.
8. Transfer max. 800 µl supernatant onto the spin filter.
9. Centrifuge for 1 min at 14,00 x g.
10. Discard spin filter.

Note: Depending on the sample a smaller amount of liquid maybe obtained after centrifugation, in this case take the actual volume.

Binding of the DNA to a Spin Filter

11. Add 200 µl **Binding Buffer** to the filtrate and mix the sample.
12. Place a new **Spin Filter** into a new yellow **2 ml Receiver Tube**.
13. Transfer max 800 µl of the sample directly onto the filter and incubate at room temperature for 1 min.
14. Centrifuge for 1 min at 14,000 x g.
15. Discard the flow-through and place the spin filter back into the same receiver tube.

Purification of the bound DNA

16. Add 550 µl **Pre-Wash Buffer** to the spin filter and centrifuge for 1 min at 14,000 x g. Discard the flow-through and place the spin filter back into the same receiver tube.
17. Add 550 µl **Wash Buffer** to the spin filter and centrifuge for 1 min at 14,000 x g. Discard the flow-through and place the spin filter back into the same receiver tube.
18. Add 550 µl **Wash Buffer** to the spin filter and centrifuge for 1 min at 14,000 x g. Discard the flow-through and place the spin filter back into the receiver tube.

Drying of the Spin Filter

19. Remove the residual ethanol by final centrifugation for 2 min at 14,000 x g.

Elution of nucleic acids from the Spin Filter

20. Place the spin filter into a clear **1.5 ml Receiver Tube**.
21. Add 100 µl of the preheated **Elution Buffer** (65 °C) directly onto the spin filter.
22. Incubate 3 min at room temperature.
23. Centrifuge for 1 min at 10,000 x g. After centrifugation discard the spin filter.

The eluted DNA is ready-to-use for PCR analyses. The DNA can be stored for up to 24 hours at 4 °C. For a storage time of more than 24 hours it should be kept at -20 °C.

Additional Information

For any inquiries or to place an order, contact your local Biotest distributor.

Domestic customers in Germany may contact Customer Service at:

Phone: +49 (0) 6103 801-496

Fax: +49 (0) 6103 801-505

mail@biotest.de

For technical inquiries please contact our Technical Support at:

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Biotest (Schweiz) AG
Biotest Italia S.r.l.
Biotest Diagnostics Corporation